



**Commercial Therapeutic/Exercise Bars with Flanges**



Specification Number:	<b>SRS-609</b>			
Product Description:	<b>Commercial Therapeutic/Exercise Bars with Flanges</b>			
S.R. Smith Part Number:	<b>TR-100A, TR-101A</b>	Scale:	NTS	
Revision:	A	Date:	07/21/03	
Written By:	T. Weil	Date:	07/21/03	
Approved By:	B. Bisceglia	Date:	07/21/03	
Page 1 of 2	Proprietary and confidential: The information contained in this specification is the sole property of S.R. Smith LLC. Any reproduction in part or whole without the written consent of S.R. Smith LLC is prohibited.			

**S.R. Smith, LLC**  
**P.O. Box 400, 1017 SW Berg Parkway**  
**Canby, Oregon 97013**  
**Tel: (888) 677-7776 (503) 266-2231**  
**Fax: (503) 266-4334**  
**www.srsmith.com**



# S.R. Smith LLC Commercial Division Product Specifications - Written

**Product Name:** Commercial Therapeutic/Exercise Bars with Flanges

**Part Number:** TR-100A, TR-101A

**General Notes:**

The **Commercial Therapeutic/Exercise Bars with Flanges** are constructed from polished stainless steel. Flanges are included to mount the bar to the pool wall. All mounting hardware is included. These bars are available in overall lengths of 48" and 88" (1219mm and 2235 mm) and special order rails in alternate lengths are also available.

**Application:**

These units are intended for use at commercial and competitive installations as exercise, training, and therapeutic rails where a 1.90" (48.3mm) diameter bar is desired. The rails can be used in most geographic areas.

**Specifications:**

- 1 Bar Assembly: These exercise bars are constructed from a 1.90" O.D. x .065 wall (48.3mm x 1.65) type 304 polished stainless steel main tube. 1.90" O.D. x .065 wall (48.3mm x 1.65) mounting posts welded to the main tube allow for a nominal clearance of 3.0" (76.2mm) between the main tube and the pool wall. Stainless steel end caps are continuously welded to the ends of the main tube to prevent corrosion.
- 2 Flange: The included flanges are made from .117" thick (2.97mm) type 304 polished stainless steel. They measure 4.42" (112.3mm) in diameter and raise .41" (10.4mm) beyond the wall. The flanges are drilled to accept three (3) included 3/8-16 UNC mounting bolts.

**Additional Product Numbers and Shipping Weight:**

Model Number	Description	Shipping Weight
TR-100A	Commercial Therapeutic/Exercise Bar with two (2) Flanges	6.0 lbs (2.72kg)
TR-101A	Commercial Therapeutic/Exercise Bar with three (3) Flanges	12.0 lbs (5.44kg)

Specification Number:	<b>SRS-609</b>		<b>S.R. Smith, LLC</b> <b>P.O. Box 400, 1017 SW Berg Parkway</b> <b>Canby, Oregon 97013</b> <b>Tel: (888) 677-7776 (503) 266-2231</b> <b>Fax: (503) 266-4334</b> <b>www.srsmith.com</b>
Product Description:	<b>Commercial Therapeutic/Exercise Bars with Flanges</b>		
S.R. Smith Part Number:	<b>TR-100A, TR-101A</b>		
Revision:	A	Date: 07/21/03	
Written By:	T. Weil	Date: 07/21/03	
Approved By:	B. Bisceglia	Date: 07/21/03	
Page 2 of 2	Proprietary and confidential: The information contained in this specification is the sole property of S.R. Smith LLC. Any reproduction in part or whole without the written consent of S.R. Smith LLC is prohibited.		